

## Aspire Overbed Table - Gunmetal Frame



- Lightweight mobile table positions overbed or chair
- Can be used for eating, drawing or other activities
- Flat top ideal for hospital, facility or home use
- Full height adjustable with a fixing knob for additional security
- C base shape fitted with four castors enables close positioning to the bedside or chair
- Powder coated mild steel frame construction in Gunmetal Grey
- Tabletop 800(W) x 430(L) mm
- Internal base 670 mm



| Code         | Display Name  | Total Height  | Colour        |
|--------------|---|---------------|---------------|
| NZ-BEA006705 | Overbed Table - Gunmetal Frame - Laminate Top - Sublime Teak  | 730 - 1095 mm | Sublime Teak  |
| NZ-BEA006706 | Overbed Table - Gunmetal Frame - Laminate Top - Milano Walnut | 730 - 1095 mm | Milano Walnut |