




VOGUE CARBON FIBRE SEAT WALKER

One of the lightest seat walkers on the market! Our Vogue Carbon Fibre Seat Walker epitomises function and style with a super lightweight carbon fibre frame for a superior walking experience.
CARBON FIBRE FRAME : Extremely lightweight to maximise portability.


COLOURS:



Small - Sapphire



Medium - Black



Tall - Emerald



LESS THAN
11lbs*



Signature Carbon Fibre Weave



1 YEAR WARRANTY

SWL
150kg
23.6st

✓ AS / NZS Approved

✓ Vogue is registered in the ARTG, AUST R ARTG 271255

- **Ergonomic Design** offers personalised height adjustment which promotes correct posture
- **Handbrakes** lock for user safety and stability
- **Simple** folding mechanism locks for easy transportation and storage
- **Bonus** Walking cane holder

- **Compact** turning circle for improved mobility in tight spaces
- **Backrest** deep seat and wide backrest for improved comfort and support
- **Storage** removable bag carries up to 5kg of essentials

CARBON FIBRE SEAT WALKER	Overall Length mm/in	Overall Width mm/in	Overall Height mm/in	Seat Depth mm/in	Seat Width mm/in	Seat Height mm/in	Folded Width mm/in	Wheel Size mm/in	Minimum Functional Weight kg/lbs	Typical Functional Weight** kg/lbs	Bag SWL kg/lbs
Small WAF705340	655 / 26	600 / 24	715 – 830 / 28 – 32	240 / 9.5	460 / 18	500 / 20	230 / 9	203 / 8	4.7 / 10.5	5.9 / 13	5 / 11
Medium WAF705350	655 / 26	615 / 24	725 – 940 / 28 – 37	240 / 9.5	460 / 18	550 / 22	230 / 9	203 / 8	4.9 / 11	6.2 / 13.7	5 / 11
Tall WAF705370	655 / 26	600 / 24	780 – 1030 / 30 – 40	240 / 9.5	460 / 18	620 / 25	230 / 9	203 / 8	5.2 / 11.5	6.3 / 14	5 / 11

* Carbon Fibre Walker - Medium | ** Includes accessories, bag and backstrap.