



WALKER RANGE

Aspire for... Comfort, stability, support



aspire WALKERS

Aspire Walkers provide assistance, support, and relief to people who want to remain actively mobile and independent — giving confidence to users of all shapes, sizes, and mobility levels.

COMMON FEATURES & BENEFITS

✓ DVA Contracted

✓ AS / NZS Approved

1 YEAR LIMITED WARRANTY



Rustproof Aluminium Frames

For durability and stability



Adjustable Handle Heights

Facilitating optimal sizing for postural positioning



Latex-free Ergonomically Shaped Grip Handles

Incredibly comfortable for users with arthritis



Ultra-Soft Hand Brakes

With effortless palm ball lock



Easy Folding Mechanism

Facilitating compact storage and transportation



Flip up Padded Seat

For accessing clip-on removable storage bags with discreet pocket compartments



Curved Padded Backrests

For seated support and comfort



360° Swivelling Front Wheels

For easy movability



Easy Wipe

For clean surfaces



Tool Free

Assembly and adjustment

CLASSIC 8 SEAT WALKER

Handle Height Range
820 – 975mm

Taller Seat Height
575mm

8" Solid Wheels
Ease over even the most robust indoor and outdoor terrain



Available colours:

■ WAF705700BL ■ WAF705700RE

✓ Classic 8 is registered in the ARTG, AUST R ARTG 224644

CLASSIC 6 SEAT WALKER

Handle Height Range
760 – 915mm

Medium Seat Height
525mm

6" Solid Wheels
Ease of mobility over every day indoor and outdoor surfaces



Available colours:

■ WAF709600BL ■ WAF709600RE

✓ Classic 6 is registered in the ARTG, AUST R ARTG 224644

DELUXE 8 SEAT WALKER

Deluxe, In-Line Rear Wheels and Brakes

Increases space around the user's stepping zone

Handle Height Range
835 – 990mm

Taller Seat Height
590mm

8" Solid Wheels
Ease over even the most robust indoor and outdoor terrain



Available colours:

■ WAF705750BL ■ WAF705750RE

✓ Deluxe is registered in the ARTG, AUST R ARTG 224650

MINI 6 SEAT WALKER

Handle Height Range
745 – 905mm

Lower Seat Height
485mm

6" Solid Wheels
Ease over every day surfaces indoors and out



Available colours:

■ WAF705800BL ■ WAF705800RE

✓ Mini 6 is registered in the ARTG, AUST R ARTG 224651

FLEX 6 SEAT WALKER

Handle Height Range
730 – 830mm

Seat Height Range
430 – 520mm

Adjustable seat and handle height

6" Solid Wheels
Ease of mobility over every day indoor and outdoor surfaces



Available colour:

■ WAF709660CH

✓ Flex 6 is registered in the ARTG, AUST R ARTG 224644

FLEX 8 SEAT WALKER

Handle Height Range
780 – 970mm

Seat Height Range
475 – 570mm

Adjustable seat and handle height

8" Solid Wheels
Ease over even the most robust indoor and outdoor terrain



Available colour:

■ WAF709670BK

✓ Flex 8 is registered in the ARTG, AUST R ARTG 224644

XL 8 SEAT WALKER

Handle Height Range
850 – 985mm

Wide Padded Seat
With 180kg safe working load for bigger-sized users

Medium Seat Height
560mm

8" Solid Wheels
Ease over even the most robust indoor and outdoor terrain



Available colour:

■ WAF750020SL

✓ XL 8 is registered in the ARTG, AUST R ARTG 224645

KNEE WALKER

Mobility Support
For an injury below the knee

Knee Pad
Adjustable height

8" Wheels
For outdoor manoeuvrability

Fold-Down Mechanism
For efficient storage



Available colour:

■ WAF712200

✓ Knee Walker is registered in the ARTG, AUST R ARTG 306095

aspire WALKERS

TRI WHEEL WALKER

Handle Height Range 825 to 925mm

Ideal for use in confined spaces

Robust steel construction

Folds for easy storage and transportation

Storage bag, basket and meal trays

8" Solid Wheels provide smooth manoeuvrability over variable surfaces



Available colour:

■ WAF705500

Tri Wheel is registered in the ARTG, AUST R ARTG 281703

MEAL TRAY WALKER

Ideal for transporting meals around the home

Clip On Plastic Trays with outer lip edge to capture spillage and prevent items from slipping off

Height Adjustable Handles with lockable handbrakes to hold the walker in place

Large 6" Wheels easier to navigate on carpets and door thresholds

360° Swivelling front wheels



Available colour:

■ DLG272061

Meal Walker is registered in the ARTG, AUST R ARTG 271256

MEAL TROLLEY

Ideal for transporting meals around the home

Clip On Plastic Trays with outer lip edge to capture spillage and prevent items from slipping off

3" Castors with two locking castors to hold the trolley in place



Available colour:

■ DLG272060

Meal Trolley is registered in the ARTG, AUST R ARTG 271256

*Please note, not all common features and benefits are applicable to these products

SPECIFICATIONS

All measurements shown are in millimetres

WALKERS	OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	SEAT DEPTH	SEAT WIDTH	SEAT HEIGHT	BACKREST HEIGHT	FOLDED WIDTH	WHEEL SIZE	UNIT WEIGHT	SWL
CLASSIC 8	820 - 975	660	630	305	355	575	815	200	8"	7kg	130kg
CLASSIC 6	760 - 915	625	615	305	355	525	770	200	6"	6kg	130kg
DELUXE 8	835 - 990	765	605	305	355	590	840	245	8"	7kg	130kg
FLEX 6	730 - 840	630	550	448	448	436 - 524	720	220	6"	6.5	136kg
FLEX 8	777 - 975	485	589	485	450	483 - 572	758	293	8"	7.5	136kg
MINI 6	745 - 905	560	615	305	355	485	690	200	6"	6kg	130kg
XL 8	850 - 985	700	735	340	495	560	850	230	8"	8kg	180kg
TRI WHEEL	825 - 925	575	715	N/A	N/A	N/A	N/A	275	8"	5.5kg	115kg

WALKERS	OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	KNEE PAD DEPTH	KNEE PAD WIDTH	KNEE PAD HEIGHT	HANDLE HEIGHT	WHEEL SIZE	UNIT WEIGHT	SWL
KNEE WALKER	840 - 970	xxx	430	360	180	480 - 560	840 - 970	8"	10kg	145kg

MEAL WALKING AIDS	OVERALL HEIGHT	OVERALL LENGTH	OVERALL WIDTH	LOWER TRAY HEIGHT	UPPER TRAY HEIGHT	TRAY WIDTH	WHEEL SIZE	UNIT WEIGHT	SWL
MEAL TRAY WALKER	770 - 920	560	585	415	770 - 920	450	6"	8kg	125kg
MEAL TROLLEY	810 - 940	505	460	355 - 480	615 - 740	450	3"	6kg	N/A

ACCESSORIES



OXYGEN BOTTLE HOLDER

Compatible with Classic/Mini/Deluxe Walkers
Attaches to A-Frame Seat Walkers
WAA686905



SEAT WALKER BASKET

Clamp On Basket To Suit Aspire Classic, Mini and Push down Seat Walkers
WAA686909



WALKING STICK HOLDER

Attaches To A-Frame Seat Walkers
Top mount clamp to securely hold the walking stick in place
WAA686924