



aspire

BATHROOM RANGE

aspire for...

Safety

Comfort

Support

Shower Chairs

Shower Stools

Over Toilet Aids

Toilet Surround Frame

Raised Toilet Seats



FEATURING
THE NEW
MAXI
RANGE

aidacare.co.nz
09 520 2962






✓ DVA Approved

✓ AS/NZS Approved

✓ Shower Chairs are registered in the ARTG, AUST # ARTG 230199/230165



Shower Chair - Wide

5 Year Warranty Corrosion

1 Year Warranty Structure


MODELS

- BTS118000 Shower Chair with Arms - Steel
- BTS118005 Shower Chair with Arms - Aluminium
- BTS118010 Shower Chair with Arms - Wide - Steel
- BTS118015 Shower Chair with Arms - Wide - Aluminium

1	Safe Working Load (S.W.L) of 200kg. Tested to AS/NZS3973.7:2009
2	Contoured plastic seat with drainage holes and windows - Provides postural support, improves drainage and reduces skin adhesion
3	Height adjustable legs - Easy push button height adjustment for changing needs
4	Additional height increments on one leg - Accommodates uneven surfaces

5	All frames are individually treated for corrosion resistance
6	Raised armrests - Assist the user transferring on and off the chair
7	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability
8	Heavy duty, reinforced frame accommodates very high SWL

SIZE	OVERALL HEIGHT	OVERALL WIDTH	OVERALL DEPTH	SEAT DEPTH	SEAT WIDTH	BETWEEN ARMS	SEAT HEIGHT	UNIT WEIGHT	S.W.L
470mm	735-910 mm	512 mm	570 mm	400 mm	440 mm	470 mm	375-550 mm	Alum: 3.1 kg Steel: 3.4 kg	200 kg
525mm (wide)	735-910 mm	570 mm	570 mm	400 mm	440 mm	525 mm	375-550 mm	Alum: 3.3 kg Steel: 3.9 kg	200 kg



✓ DVA Approved

✓ AS/NZS Approved

✓ MAXI is registered in the ARTG, AUST # ARTG 304118

5 Year Warranty Corrosion

1 Year Warranty Structure

MODEL

BTS118020 | MAXI Shower Chair

1	Safe Working Load (S.W.L) of 310kg.
2	Super heavy duty width adjustable frame to 775mm
3	Soft padding on seat and arms for postural support and comfort
4	Reinforced aluminium frame accommodates bariatric users and is corrosion resistant

5	Height adjustable legs with easy circlip adjustment points for changing needs
6	Raised armrests - Assist the user transferring on and off the chair
7	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability


OVERALL HEIGHT	OVERALL WIDTH	OVERALL DEPTH	SEAT DEPTH	SEAT WIDTH	BETWEEN ARMS	SEAT HEIGHT	S.W.L
670-1005 mm	645-845 mm	495 mm	455 mm	575-775 mm	575-775 mm	455-575 mm	310 kg



- DVA Approved
- AS/NZS Approved
- Shower Stools are registered in the ARTG, AUST R ARTG 230166/230167



BTS119025



7

MODELS

- BTS119000 Shower Stool - Steel
- BTS119005 Shower Stool - Aluminium
- BTS119020 Shower Stool with Padded Seat (PS) - Steel
- BTS119025 Shower Stool with Padded Seat & Back (PS & B) - Steel



1	Safe Working Load (S.W.L) of 175kg. Tested to AS/NZS3973.7:2009
2	Clip on seat or padded back attachment – For easy assembly and cleaning
3	Height adjustable legs – Easy push button height adjustment for changing needs
4	All frames are individually treated for corrosion resistance

5	Seat with anterior tilt – Facilitates standing from sitting and provides greater comfort for stiff or painful hips
6	Raised armrests – Assist the user transferring on and off the stool
7	Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
8	Additional height increments on one leg – Accommodates uneven surfaces

OVERALL HEIGHT	SEAT HEIGHT	OVERALL WIDTH	OVERALL DEPTH	SEAT DEPTH	SEAT WIDTH	BETWEEN ARMS	UNIT WEIGHT	S.W.L
670-855 mm	470-655 mm	545 mm	400 mm	310 mm	375 mm	495 mm	Alum: 3.1 kg Steel: 3.8 kg	175 kg



- DVA Approved
- AS/NZS Approved
- Shower Chairs are registered in the ARTG, AUST R ARTG 230166/230167



2



7

MODELS

- BTS119010 Shower Stool - Wide - Steel
- BTS119015 Shower Stool - Wide - Aluminium



1	Safe Working Load (S.W.L) of 175kg. Tested to AS/NZS3973.7:2009
2	Clip on seat attachment – For easy assembly and cleaning
3	Height adjustable legs – Easy push button height adjustment for changing needs
4	All frames are individually treated for corrosion resistance
5	Seat with anterior tilt – Facilitates standing from sitting and provides greater comfort for stiff or painful hips

6	Raised armrests - Assist the user transferring on and off the stool
7	Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
8	Additional height increments on one leg – Accommodates uneven surfaces
9	Comfortable wide seat with drainage holes – For increased support and hygiene
10	Reinforced wide frame – Provides stability for larger users

OVERALL HEIGHT	SEAT HEIGHT	OVERALL WIDTH	OVERALL DEPTH	SEAT DEPTH	SEAT WIDTH	BETWEEN ARMS	UNIT WEIGHT	S.W.L
575-765 mm	425-615 mm	570 mm	500 mm	365 mm	420 mm	525 mm	Alum: 3.1 kg Steel: 3.8 kg	175 kg



- ✓ DVA Approved
- ✓ AS/NZS Approved
- ✓ MAXI is registered in the ARTG, AUST # ARTG 206118




5 Year Warranty
Corrosion

1 Year Warranty
Structure

MODEL
BTS119030 | MAXI Shower Stool

1	Safe Working Load (S.W.L) of 310kg.	5	Height adjustable legs with easy circlip adjustment points for changing needs
2	Super heavy duty width adjustable frame to 775mm	6	Raised armrests – Assist the user transferring on and off the stool
3	Soft padding on seat and arms for postural support and comfort	7	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability
4	Reinforced aluminium frame accommodates bariatric users and is corrosion resistant		



- ✓ DVA Approved
- ✓ AS/NZS Approved
- ✓ Toilet Surround Frame is registered in the ARTG, AUST # ARTG 230170




5 Year Warranty
Corrosion

1 Year Warranty
Structure

MODEL
BTT155100 | Toilet Surround - Aluminium

1	Safe Working Load (S.W.L) of 175kg. Tested to AS/NZS3973.7:2009	4	Durable non-slip feet provide excellent grip and stability in wet conditions and a stainless steel washer enhances durability
2	Height adjustable frame – Easy push button height adjustment for changing needs	5	Additional height increments on one leg – Accommodates uneven surfaces
3	Raised armrests – Assist the user transferring on and off the toilet		





- ✓ DVA Approved
- ✓ AS/NZS Approved
- ✓ Over Toilet Aids are registered in the ARTG, AUST # ARTG 230169/230170

Over toilet Aid - Wide

MODELS

BTT145000	Over Toilet Aid - Steel
BTT145005	Over Toilet Aid - Aluminium
BTT145010	Over Toilet Aid - Wide - Steel
BTT145015	Over Toilet Aid - Wide - Aluminium



ACCESSORY

Aspire Dementia Care Toilet Seat
BTT145102



1	Safe Working Load (S.W.L) of 175kg. Tested to AS/NZS3973.7:2009
2	Raised armrests – Assist the user transferring on and off the over toilet aid
3	Height adjustable legs – Easy push button height adjustment for changing needs
4	Comfortable, contoured seat with clip on attachment – For easy assembly and cleaning. Padded seat available, sold separately (BTT145110)

5	All frames are individually treated for corrosion resistance
6	Removable, moulded splash guard provided – Easily removed for transport and cleaning
7	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability
8	Additional height increments on one leg – Accommodates uneven surfaces

SIZE	OVERALL HEIGHT	OVERALL WIDTH	OVERALL DEPTH	SEAT DEPTH	SEAT WIDTH	BETWEEN ARMS	SEAT HEIGHT	UNIT WEIGHT	S.W.L
460 mm	700-880 mm	512 mm	625 mm	435 mm	370 mm	460 mm	450-630 mm	Alum: 3.4 kg Steel: 4.2 kg	175 kg
525 mm (Wide)	700-880 mm	570 mm	625 mm	435 mm	370 mm	525 mm	450-630 mm	Alum: 3.4 kg Steel: 4.2 kg	175 kg



- ✓ DVA Approved
- ✓ AS/NZS Approved
- ✓ MAXI is registered in the ARTG, AUST # ARTG 306118

MODEL

BTT145020 | MAXI Over Toilet Aid

1	Safe Working Load (S.W.L) of 310kg.
2	Super heavy duty width adjustable frame to 775mm
3	Flip up padded seat for ease of cleaning
4	Reinforced aluminium frame accommodates bariatric users and is corrosion resistant

5	Height adjustable legs with easy circlip adjustment points for changing needs
6	Raised armrests – Assist the user transferring on and off the over toilet aid
7	Durable non-slip feet provide excellent grip and stability in wet conditions. Stainless steel washer enhances durability



OVERALL HEIGHT	OVERALL WIDTH	OVERALL DEPTH	SEAT DEPTH	SEAT WIDTH	BETWEEN ARMS	SEAT HEIGHT	S.W.L
670-790 mm	645-845 mm	495 mm	455 mm	575-775 mm	575-775 mm	455-575 mm	310 kg



- ✓ DVA Approved
- ✓ AS/NZS Approved
- ✓ Toilet Seat Raisers are registered in the ARTG, AU21 & ARTG 242921

MODELS

BTT149800	Toilet Seat Raiser - 50mm - No Lid
BTT150200	Toilet Seat Raiser - 100mm - No Lid
BTT150300	Toilet Seat Raiser - 150mm - No Lid
BTT149900	Toilet Seat Raiser - 50mm - With Lid
BTT150100	Toilet Seat Raiser - 100mm - With Lid
BTT150400	Toilet Seat Raiser - 150mm - With Lid

5 Year Warranty
Corrosion


1 Year Warranty
Structure

- | | | | |
|---|--|---|---|
| 1 | Three heights available – helps facilitate sit-to-stand transfers for safe and comfortable toileting | 4 | Fully sealed structure – Easy to clean for increased infection control |
| 2 | Contoured plastic seat – Comfortable and durable | 5 | Generous cut outs at the front and rear – Facilitates personal care and hygiene |
| 3 | Available with or without lid | 6 | Two adjustable mounting brackets – Provide an adjustable, secure fit |


OVERALL HEIGHT	UNIT WEIGHT	S.W.L
2" (50mm)	1.4 kg	200 kg
4" (100mm)		
6" (150mm)		



SHOWER COMMODES



BATH SPONGES



BATH LIFTS



BATH/SHOWER TRANSFER BENCHES



BATH STEPS



GRAB RAILS



HAND HELD SHOWERS



SWIVEL BATHERS



BATH BOARDS

CRUTCHES, FRAMES & STICKS



SCOOTERS



OVERBED TABLES



BED STICKS & POLES



RAMPS



AIR MATTRESSES & FOAM MATTRESSES



COMMODOES, URINALS & BED PANS



LIFTERS & SLINGS



BEDS & BED POLES



DAILY LIVING AIDS



SHOWER COMMODOES



BATH BOARDS & BENCHES



SEAT WALKERS



POWER WHEELCHAIRS



PRESSURE CUSHIONS



LIFT & RECLINE CHAIRS



Distributed by:

