




# DAY CHAIRS, STOOLS & COMMODOES

Aspire for... Comfort, stability, support




# HIGH BACK CLASSIC DAY CHAIR

A high back orthopaedic chair that offers adjustability, comfort and support.



Height Adjustable



Postural Support




## Contoured High Backrest

With a slight recline angle provides superior postural support and opens hip positioning for comfort



Also Available in Slate Vinyl



## MODELS

- CHP208050 High Back Classic Day Chair – Champagne Vinyl
- CHP208055 High Back Classic Day Chair – Slate Vinyl

# WATERFALL DAY CHAIR

A high back orthopaedic chair with contoured and adjustable backrest cushions that offer improved comfort and postural support.



Ink Vinyl




Latte Vinyl



Oatmeal Fabric



Mocha Fabric



## Waterfall

Tri-Pillow Backrest Provides comfort and support

## Removable Cushions

Seat and Back Cushion overlays can be easily detached for cleaning

## Upholstery Range

Ink and Latte vinyl plus Oatmeal and Mocha fabric options




Height Adjustable



## MODELS

- CHP208070 Waterfall Day Chair – Ink Vinyl
- CHP208071 Waterfall Day Chair – Latte Vinyl
- CHP208075 Waterfall Day Chair – Oatmeal Fabric
- CHP208076 Waterfall Day Chair – Mocha Fabric

## FEATURES & BENEFITS SHARED BY THE HIGHBACK CLASSIC AND WATERFALL DAY CHAIRS



### Height Adjustable Legs

With easy push button adjustment for comfortable positioning and changing needs



### Sturdy

**Lightweight Design**  
For a generous weight capacity, yet lightweight for easy manoeuvrability



### Non-Slip Tips

On chair legs provide stability whilst sitting or transferring



### Robust Powder Coated Steel Frames

Are durable and easy to clean



### Softly Padded Armrests, Seat & Back

Distribute pressure whilst simultaneously offering support surfaces for prolonged sitting comfort




### Transport Legs With Wheels

Available as an accessory to facilitate easier manoeuvring of the chairs (without a patient in the chair)

See back cover for details



### Easy Clean Vinyl Upholstery



✓ DVA CONTRACTED

✓ AS / NZS Approved

✓ Waterfall Day Chair is registered in the ARTG, AUST R ARTG 230087

✓ High Back Classic Day Chair is registered in the ARTG, AUST R ARTG 230087

Seat Depth (mm)	Overall Depth (mm)	Seat Width (mm)	Overall Width (mm)	Seat Height (mm)	Backrest Height (mm)	Overall Height (mm)	Unit Weight (kg)	SWL (kg)
490	600	490	645	425-610	720	1160-1312	16	160

Seat Depth (mm)	Overall Depth (mm)	Seat Width (mm)	Overall Width (mm)	Seat Height (mm)	Backrest Height (mm)	Overall Height (mm)	Unit Weight (kg)	SWL (kg)
490	670	535	720	420-580	685	1045-1200	16.5	200

# ADJUSTABLE & PRESSURE REDUCING DAY CHAIRS

A truly versatile height, width and depth adjustable high back orthopaedic chair with built in lumbar and lateral support and knee bolster. Available in pressure reducing variety.

## ADJUSTABLE DAY CHAIR

**Head Cushion**  
Height adjustable and removable head cushion for extra comfort

**Contoured Seat and Backrest Padding**  
Designed for pressure redistribution and user comfort

**Sturdy design**  
For a generous weight capacity of 160kg, yet lightweight for easy manoeuvrability



Also Available in Latte Vinyl



510 - 610 mm

## PRESSURE REDUCING ADJUSTABLE DAY CHAIR

**Pressure Redistribution**  
Carrflex™ PU 4-Way Stretch Fabric for superior chemical resistance, retains flexibility and breathability

**High quality Australian memory foam** provides a comfortable and supportive surface, improving user comfort and pressure redistribution

**Maximised Immersion and Envelopment** to effectively redistribute pressure, improves blood flow and promotes healing



Available in 4-way stretch fabric



## MODELS

- CHP208025 Adjustable Day Chair – Latte Vinyl
- CHP208030 Adjustable Day Chair – Ink Vinyl

- CHP208035 Adjustable Day Chair - Pressure Reducing Blue Stretch Material

# MAXI ADJUSTABLE DAY CHAIR

Bariatric Day Chair with 310kg weight capacity is fully height, width and depth adjustable to suit an array of large user shapes and sizes

**Seat with Anterior Tilt**  
Facilitates standing from sitting and provides greater comfort for stiff or painful hips

**Reinforced Design and Weight Capacity**  
Caters to a MAXI weight capacity of 310kg especially for bariatric users



600-810mm



Push Handle & Transport Wheels



Adjustable Seat Depth & Backrest Height




Height & Width Adjustable Frame

## MODELS

- CHP208010 Maxi Adjustable Day Chair – Burgundy Vinyl

## FEATURES & BENEFITS SHARED BY THE ADJUSTABLE AND MAXI ADJUSTABLE DAY CHAIRS




**Height Adjustable Legs**  
With easy push button adjustment for comfortable positioning and changing needs



**Tool Free Adjustment Points**  
Allow the chair to be customised to suit long term use or readily adapted to suit multiple users



**Adjustable Seat Depth, Width and Height**  
Allows customised configuration to a user's needs and circumstances



**Unique Width Adjustable Frames**  
Cater to the full spectrum of users shapes and sizes



**Softly Padded Armrests**  
Distribute pressure whilst simultaneously offering support surfaces for prolonged sitting comfort



**Push Handle and Transport Wheels**  
Allow a carer to tilt the chairs backwards onto the rear leg wheels in order to manoeuvre it about (without a patient in the chair)

See back cover for details



**Easy Clean Vinyl Upholstery**

1 YEAR WARRANTY

✓ DVA CONTRACTED

✓ AS / NZS Approved

✓ Maxi Day Chair is registered in the ARTG, AUST R ARTG 230087

✓ Adjustable Day Chair is registered in the ARTG, AUST R ARTG 230087

Seat Depth (mm)	Overall Depth (mm)	Seat Width (mm)	Overall Width (mm)	Seat Height (mm)	Backrest Height (mm)	Overall Height (mm)	Unit Weight (kg)	SWL (kg)
475-590	760-610	510-610	575-675	410-530	740-790	1070-1220	26	160

Seat Depth (mm)	Overall Depth (mm)	Seat Width (mm)	Overall Width (mm)	Seat Height (mm)	Backrest Height (mm)	Overall Height (mm)	Unit Weight (kg)	SWL (kg)
390-530	850-990	600-810	600-810	370-520	1120-1340	1440	45	310



## LOWBACK DAY CHAIR

A versatile, low back orthopaedic chair



**Contoured and Reclined Backrest**  
Contour and slight recline of backrest angle helps maintain posture while supporting an open and comfortable hip position



Also Available in Champagne Vinyl

### MODELS

- CHP208000 Low Back Classic Day Chair – Champagne Vinyl
- CHP208005 Low Back Classic Day Chair – Slate Vinyl

## PERCHING STOOL

A height adjustable stool that provides relief from fatigue in the kitchen or living area



**Raised Armrests**  
Assist with sitting and standing

**Seat with Anterior Tilt**  
Facilitates easier standing and sitting and provides greater comfort for stiff or painful hips



### MODELS

- CHP213500 Perching Stool – Champagne Vinyl

## CLASSIC BEDSIDE COMMUNE

Allows for toileting in the bed room for users unable to travel between the bedroom and bathroom



**Padded Armrests**  
Are raised to provide an easily reachable support point for lowering onto the seat

**Lightweight**  
Bedside commode with padded seat, back and arms



Also Available in Slate Vinyl

### MODELS

- BEC045005 Classic Bedside Commode – Champagne Vinyl
- BEC045006 Classic Bedside Commode – Slate Vinyl

## DELUXE BEDSIDE COMMUNE

Allows for toileting in the bed room with increased comfort for users unable to travel between the bedroom and bathroom



**Contoured and Reclined Backrest**  
Contour and slight recline of backrest angle helps maintain posture while supporting an open and comfortable hip position

**Softly Padded Armrests, Seat and Back**  
Offer support surfaces for prolonged sitting comfort




Also Available in Champagne Vinyl

### MODELS

- BEC045115 Deluxe Bedside Commode – Champagne Vinyl
- BEC045116 Deluxe Bedside Commode – Slate Vinyl

### FEATURES & BENEFITS SHARED BY THE LOWBACK DAY CHAIR AND PERCHING STOOL



**Height Adjustable Legs**  
Easy push button adjustment for comfortable positioning



**Soft Padding**  
Distributes pressure and offers support for prolonged sitting comfort



**Non-Slip Tips**  
Provides stability whilst sitting or transferring



**Robust Powder Coated Steel Frames**  
Are durable and easy to clean



**Easy Clean**  
Vinyl upholstery



**Lightweight**  
For easy manoeuvrability

**1 YEAR WARRANTY**

**DVA CONTRACTED**


Classic Bedside Commode is registered in the ARTC, AUST R ARTC 230171

**AS / NZS Approved**

Deluxe Bedside Commode is registered in the ARTC, AUST R ARTC 230171

	Seat Depth (mm)	Overall Depth (mm)	Seat Width (mm)	Overall Width (mm)	Seat Height (mm)	Backrest Height (mm)	Overall Height (mm)	Unit Weight (kg)	SWL (kg)
Low Back Chair	440	500	460	580	450-605	440	865-1015	9	160
Perching Stool	310	470	375	550	640-690	440	940-1090	6.8	175

### FEATURES & BENEFITS SHARED BY THE CLASSIC AND DELUXE BEDSIDE COMMUNES



**Height Adjustable Legs**  
Easy push button adjustment for comfortable positioning



**Seat Access**  
Folds out of the way to access the toilet seat and commode bowl



**Three-Sided Shroud**  
Discreetly conceals bowl



**Lightweight**  
For easy manoeuvrability



**Removable Bowl**  
With lid and carry handle



**Easy Clean**  
Vinyl upholstery

**1 YEAR WARRANTY**

**DVA CONTRACTED**

Classic Bedside Commode is registered in the ARTC, AUST R ARTC 230171

**AS / NZS Approved**

Deluxe Bedside Commode is registered in the ARTC, AUST R ARTC 230171

	Seat Depth (mm)	Overall Depth (mm)	Seat Width (mm)	Overall Width (mm)	Seat to Floor Height (mm)	Toilet Seat to Floor Height (mm)	Backrest Height (mm)	Overall Height (mm)	Unit Weight (kg)	SWL (kg)
Classic	440	570	450	570	500-625	475-600	455	915-1040	5.5	125
Deluxe	450	520	460	600	480-630	435-585	425	865-1015	12	160

# aspire DAY CHAIRS AND STOOLS ACCESSORIES

## TRANSPORT LEGS WITH WHEELS

✓ Transport Legs with Wheels is registered in the ARTG, AUST R ARTG 230087

### FEATURES & BENEFITS

- Assists carers to safely and easily tilt and manoeuvre an empty chair around a room
- Wheeled legs insert directly into the frame to replace regular rear legs
- Transport wheels come as standard with the Adjustable day chair and Maxi Adjustable Day Chair and can be fitted as options on the Waterfall and High Back Classic day chairs



Transport Wheels  
CHP208090



### COMPATIBLE CHAIRS

- CHP208050 | High Back Classic Day Chair – Champagne Vinyl
- CHP208055 | High Back Classic Day Chair – Slate Vinyl
- CHP208070 | Waterfall Day Chair – Ink Vinyl

- CHP208071 | Waterfall Day Chair – Latte Vinyl
- CHP208075 | Waterfall Day Chair – Oatmeal Fabric
- CHP208076 | Waterfall Day Chair – Mocha Fabric

## PADDED LEGREST

✓ Aspire Padded Legrest is registered in the ARTG, AUST R ARTG 230087



Champagne Vinyl



Slate Vinyl


- CHP208080 | Aspire Padded Legrest – Champagne Vinyl
- CHP208085 | Aspire Padded Legrest – Slate Vinyl

### FEATURES & BENEFITS



#### Comfortable Resting Solution

Legs can be raised when sitting in a day chair to provide extra comfort




#### Height Adjustable

Easy push button adjustment for comfortable positioning and changing needs



#### Robust Powder Coated Steel Frames

Are durable and easy to clean



#### Padded Surface

Generous support of legs and feet