



# EVOKE 2 WHEELCHAIR

The Aspire EVOKE 2 and EVOKE 2 HD are highly adjustable, lightweight, aluminium wheelchairs designed to enable a personalised fit for every user. A comprehensive range of compatible accessories is available for increased independence and enhanced clinical outcomes.



CRASH TEST  
**APPROVED**  
AS/NZS 3696.19

**1** YEAR WARRANTY  
Components

**2** YEARS WARRANTY  
Frame



**SILVER**  
FRAME



**BLACK**  
UPHOLSTERY

**EVOKE 2**



**EVOKE 2 HD**



**Swing-In / Out, Height Adjustable & Removable Foot Rests**

This facilitates user transfers and allows optimal positioning for the changing needs of the user

**Strong, Nylon Upholstery with a Comfortable Cushion & Padded Backrest**

For increased comfort, a padded overlay is available as an accessory

**Depth Adjustable Seating with Rear Angle Adjustability**

For more accurate positioning, customisability and increased user comfort

**Quick Release 24" Solid Rear Tyres**

The wheels are easily removed to reduce the weight of the chair for transport and storage



\*Scan here for more information on the Evoke 2 Wheelchair Attendant Brake options

Pictured with optional transit conversion kit



**ADJUSTMENTS**



**Seat & Backrest Height**



**Seat Depth**



**Backrest Angle**

**OPTIONAL ACCESSORIES**



**Contoured Padded Overlay**



**Attendant Brake Kit & Conversion Kit\***

**EVOKE 2 WHEELCHAIR**

	Seat Width mm	Seat Depth mm	Seat Height mm	Overall Width mm	Backrest Height mm	Armrest Height mm	Min. Chair Weight kg	SWL kg
MWS449500 Standard 350mm	350	300 – 450	450 – 500	560	430 – 490	210 – 230	14.8	150
MWS449510 Standard 400mm	400	300 – 450	450 – 500	610	430 – 490	210 – 230	14.9	150
MWS449520 Standard 450mm	450	300 – 450	450 – 500	660	430 – 490	210 – 230	15	150
MWS449530 Standard 500mm	500	400 – 550	450 – 500	710	430 – 490	210 – 230	15.1	150
MWS449540 HD 570mm	570	400 – 550	450 – 500	780	430 – 490	210 – 230	16	190
MWS449550 HD 610mm	610	400 – 550	450 – 500	820	430 – 490	210 – 230	16.2	190

**ACCESSORIES**

WCA848592 Attendant Brake Kit (No Wheels)

WCA848594 Transit Conversion Kit (Incl. Brake and Wheels)

\*Range of motion dimensions are accurate to a +/-10mm tolerance

ASFY123.2.221T24