aspire GO TURNER

The Aspire Go Turner is a proven sit-to-sit transfer aid for quick and easy manual handling. This compact folding patient mover has a very narrow footprint, making it the perfect solution for transferring patients in tight spaces. The pull apart and folding design reduces storage space requirements and helps to transport the product if needed.



BELT COMPATIBILITY

The Go Turner features four different sized belts for safety and comfort during transportation.

Code	Description	Hip Circumference	S.W.L
LSS391132	Aspire Go Turner Support Belt - Small	78 - 98cm	150kg
LSS391133	Aspire Go Turner Support Belt - Medium	93 - 118cm	150kg
LSS391134	Aspire Go Turner Support Belt - Large	113 - 138cm	150kg
LSS391135	Aspire Go Turner Support Belt - XLarge	133 - 158cm	150kg

*Range of motion dimensions are accurate to a +/-10mm tolerance