



WALKER RANGE

Premium Mobility Solutions




VOGUE CARBON FIBRE SEAT WALKER

One of the lightest seat walkers on the market! Our Vogue Carbon Fibre Seat Walker epitomises function and style with a super lightweight carbon fibre frame for a superior walking experience.

CARBON FIBRE FRAME: Extremely lightweight to maximise portability

COLOURS:



Signature Carbon Fibre Weave

LESS THAN 5KG*



MODELS


- WAF705340 Vogue Carbon Fibre - Small - Sapphire
- WAF705350 Vogue Carbon Fibre - Medium - Black
- WAF705370 Vogue Carbon Fibre - Tall - Emerald

VOGUE LIGHTWEIGHT 2 SEAT WALKER

This modern seat walker has a clean and crisp frame design, making it stand out when going for a walk or browsing the shops!

LIGHTWEIGHT ALUMINIUM BODY: Improves manoeuvrability and portability.

COLOURS:



MODELS

- WAF705310RE Vogue Lightweight 2 - Red
- WAF705310CH Vogue Lightweight 2 - Champagne
- WAF705310LE Vogue Lightweight 2 - Leopard
- WAF705310FL Vogue Lightweight 2 - Flower
- WAF705310ZB Vogue Lightweight 2 - Zebra




Discreet Cable Management

JUST UNDER 6KG




FEATURES & BENEFITS SHARED BY THE VOGUE CARBON FIBRE ...



Ergonomic Design
Offers personalised height adjustment which promotes correct posture



Compact
Turning circle for improved mobility in tight spaces




Simple
Folding mechanism locks for easy transportation and storage




Handbrakes
Lock for user safety and stability

... AND LIGHTWEIGHT ALUMINIUM SEAT WALKERS



Storage
Removable bag carries up to 5kg of essentials



Backrest
Deep seat and wide backrest for improved comfort and support



Bonus
Walking cane holder



AS / NZS Approved

Vogue is registered in the ARTG, AUST R ARTG 271255

	Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL
Small	655mm	600mm	715-830mm	240mm	460mm	500mm	230mm	203mm	4.7kg	5.9kg	5kg
	26"	23.6"	28"-32"	9.4"	18"	20"	9"	8"	10.4lbs	13lbs	11lbs
Medium	655mm	615mm	725-940mm	240mm	460mm	550mm	230mm	203mm	4.95kg	6.2kg	5kg
	26"	24"	28"-37"	9.4"	18"	22"	9"	8"	11lbs	13.7lbs	11lbs
Tall	655mm	600mm	780-1030mm	240mm	460mm	620mm	230mm	203mm	5.2kg	6.35kg	5kg
	26"	23.6"	30"-40"	9.4"	18"	25"	9"	8"	11.5lbs	14lbs	11lbs

	Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL
Small	655mm	620mm	780-920mm	240mm	465mm	540mm	220mm	203mm	5.4kg	6.45kg	5kg
	26"	24.5"	31"-36"	9.5"	18"	21"	8.6"	8"	11.7lbs	14.2lbs	11lbs

* Carbon Fibre Walker - Medium | ** Includes accessories, bag and backstrap.


VOGUE ADVENTURE SEAT WALKER

Rugged and durable, the Adventure Seat Walker is the perfect all terrain walker. Super sized front wheels can manage any surface, perfect for the adventurous user!

COLOURS:



Black



Wider Seat supports up to 200kg



Large Multi-Terrain Wheels



Folds to 245mm



MODELS

WAF705450 | Vogue Adventure - Black

FEATURES & BENEFITS



Ergonomic Design
Offers personalised height adjustment which promotes correct posture




Compact
Turning circle for improved mobility in tight spaces



Handbrakes
Lock for user safety and stability



Storage
Removable bag carries up to 5kg of essentials



Simple
Folding mechanism locks for easy transportation and storage



Bonus
Walking cane holder



1 YEAR WARRANTY

✓ AS / NZS Approved

✓ Vogue is registered in the ARTC, AUST R ARTC 271255

Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size (Front / Rear)	Minimum Functional Weight	Typical Functional Weight**	Bag SWL
770mm	695mm	900-1035mm	230mm	540mm	530mm	245mm	279mm / 225mm	8.8kg	10kg	5kg
30.3"	27.4"	35.4"-40.7"	9.1"	21.3"	20.9"	9.6"	11" / 8.8"	19.4lbs	22lbs	11lbs

VOGUE FOREARM WALKER

For those who need a little more support, the Forearm Walker is for you. Soft and contoured forearm pads ensure maximum comfort and a smooth ride.

LIGHTWEIGHT ALUMINIUM BODY: Improves manoeuvrability and portability.

COLOURS:




White / Black



Loosen. Raise/Lower. Tighten.



Loosen. Extend/Shorten. Tighten.



Padded Forearm Supports



Folds to 225mm




MODELS

WAF705400 | Vogue Forearm - White/Black

FEATURES & BENEFITS



Ergonomic Design
Offers personalised height adjustment which promotes correct posture



Simple
Folding mechanism locks for easy transportation and storage



Handbrakes
Lock for user safety and stability



Storage
Removable bag carries up to 5kg of essentials



Bonus
Walking cane holder



1 YEAR WARRANTY

✓ AS / NZS Approved

✓ Vogue is registered in the ARTC, AUST R ARTC 271255

Overall Length	Overall Width	Overall Height	Seat Depth	Seat Width	Seat Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL
655mm	608mm	1020-1270	240mm	441mm	620mm	225mm	203mm	9.1kg	10kg	5kg
26"	23.9"	40.2"-50"	9.4"	17.4"	24.4"	8.8"	8"	20.1lbs	22lbs	11lbs

** Includes accessories, bag and backstrap.

VOGUE INDOOR WALKER

Ideal for use within the home, the Indoor Walker is perfect for portable dining and room to room mobility support.


COLOURS:



Red / Black



Champagne / Black



Folds to 270mm

MODELS

- WAF705100RE | Vogue Indoor - Red/Black
- WAF705100CH | Vogue Indoor - Champagne/Black

FEATURES & BENEFITS



Ergonomic Design
Offers personalised height adjustment which promotes correct posture



Compact
Turning circle for improved mobility in tight spaces



Handbrakes
Lock for user safety and stability



Storage
Removable bag carries up to 5kg of essentials



Simple
Folding mechanism locks for easy transportation and storage



Bonus
Walking cane holder



Portable Dining
Safely collect and transport meals or hot drinks



Turn Around.
Slide Over Knee.
Lock & Enjoy.



1 YEAR WARRANTY

AS / NZS Approved

Vogue is registered in the ARTG, AUST R ARTG 271255

Overall Length	Overall Width	Overall Height	Folded Width	Wheel Size	Minimum Functional Weight	Typical Functional Weight**	Bag SWL	Tray SWL
650mm 25.6"	540mm 21.3"	860-970mm 33.9"-38.2"	270mm 10.7"	152mm 6"	4.4kg 9.7lbs	5.1kg 11.2lbs	5kg 11lbs	5kg 11lbs

** Includes accessories, bag and tray.

VOGUE WALKER ACCESSORIES



Oxygen Bottle Holder Bag*
WAA695600



Night Light - LED Lamp
WAA695630



IV Pole - Height Adjustable
WAA695651



Single Hand Brake System* - Compatible with Aspire Vogue Carbon Walkers
WAA695640

*Single Hand Brake System and Oxygen Bottle Holder Bag are not compatible with the Vogue Indoor Walker





MOBILITY SCOOTERS

Wide range of Aspire Scooters for all user needs.

Shown Aspire HS268 Manual Folding Mobility Scooter SCT671135



COMFIMOTION COLLECTION

Range of adjustable beds and accessories for both clinical needs and comfort at home.

Shown Aspire ComfiMotion Activ Care Bed in Portebello Grey - Long Single - BEB046740



LIFT RECLINE CHAIRS

Luxury Aspire recliner range has four, dual and single motor functions.

Shown Aspire PostureFit Configurable Lift Recline Chair
Sizes Available Medium & Tall - Medium/Large



OVERBED TABLES

Lightweight mobile table positions easily over the lifestyle community beds.

Shown Over Bed Table with Gunmetal Frame & Sublime Teak Laminate Top BEA006705